


## 4-(4-Chlorophenyl)piperidin-4-ol

## Chemical Properties

CAS No. :	39512-49-7
Formula:	C <sub>11</sub> H <sub>14</sub> ClNO
Molecular Weight:	211.688
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-(4-Chlorophenyl)piperidin-4-ol is a useful organic compound for research related to life sciences. The catalog number is Fr13199 and the CAS number is 39512-49-7.
Targets(IC50)	Others, iGluR

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7239 mL	23.6194 mL	47.2389 mL
5 mM	0.9448 mL	4.7239 mL	9.4478 mL
10 mM	0.4724 mL	2.3619 mL	4.7239 mL
50 mM	0.0945 mL	0.4724 mL	0.9448 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481