


1,3-Dihydro-2H-isoindol-2-ylacetic acid hydrochloride

Chemical Properties

CAS No. : 1215322-31-8
 Formula: C₁₀H₁₂ClNO₂
 Molecular Weight: 213.661
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,3-Dihydro-2H-isoindol-2-ylacetic acid hydrochloride, with CAS No. 1215322-31-8, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 1,3-Dihydro-2H-isoindol-2-ylacetic acid hydrochloride provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6803 mL	23.4017 mL	46.8033 mL
5 mM	0.9361 mL	4.6803 mL	9.3607 mL
10 mM	0.468 mL	2.3402 mL	4.6803 mL
50 mM	0.0936 mL	0.468 mL	0.9361 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481