


Tanshindiol B

Chemical Properties

CAS No. :	97465-70-8
Formula:	C ₁₈ H ₁₆ O ₅
Molecular Weight:	312.32
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tanshindiol B possesses a unique anti-cancer activity whose mechanism involves the inhibition of EZH2 activity and would provide chemically valuable information for designing a new class of potent EZH2 inhibitors. It also possesses the anti-angiogenic activity, so it can be a promising candidate for angiogenesis inhibitors.
Targets(IC50)	Reductase, Histone Methyltransferase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2018 mL	16.0092 mL	32.0184 mL
5 mM	0.6404 mL	3.2018 mL	6.4037 mL
10 mM	0.3202 mL	1.6009 mL	3.2018 mL
50 mM	0.064 mL	0.3202 mL	0.6404 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481