


3-Methyluridine

Chemical Properties

CAS No. :	2140-69-4
Formula:	C10H14N2O6
Molecular Weight:	258.23
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Methyluridine is a modified RNA nucleoside of cellular RNA and commonly found in the form of RNA modifications.
Targets(IC50)	Endogenous Metabolite

Solubility Information

Solubility	DMSO: 160 mg/mL (619.6 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 4 mg/mL (15.49 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8725 mL	19.3626 mL	38.7252 mL
5 mM	0.7745 mL	3.8725 mL	7.745 mL
10 mM	0.3873 mL	1.9363 mL	3.8725 mL
50 mM	0.0775 mL	0.3873 mL	0.7745 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481