


## Acedapsone

## Chemical Properties

CAS No. :	77-46-3
Formula:	C <sub>16</sub> H <sub>16</sub> N <sub>2</sub> O <sub>4</sub> S
Molecular Weight:	332.37
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Acedapsone has antimalarial and antimicrobial action, but is mainly used as a depot leprostatic agent.
Targets(IC50)	Others,Histone Methyltransferase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0087 mL	15.0435 mL	30.087 mL
5 mM	0.6017 mL	3.0087 mL	6.0174 mL
10 mM	0.3009 mL	1.5043 mL	3.0087 mL
50 mM	0.0602 mL	0.3009 mL	0.6017 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Olivares Sabournin R, González Segredo AB. [Bacteriologic study in patients with multibacillary leprosy treated with rifampicin and acedapsone]. Rev Cubana Med Trop. 1989 May-Aug;41(2):307-12. Spanish. PubMed PMID: 2486227.

George J, Balakrishnan S. Blood dapsone levels in leprosy patients treated with acedapsone. Indian J Lepr. 1986 Jul-Sep;58(3):401-6. PubMed PMID: 3794408.

Neelan PN, Sirumban P, Sivaprasad N. Limited duration acedapsone prophylaxis in leprosy. Indian J Lepr. 1986 Apr-Jun;58(2):251-6. PubMed PMID: 3543161.

Russell DA, Worth RM, Jano B, Fasal P, Shepard CC. Acedapsone in the prevention of leprosy: field trial in three high prevalence villages in Micronesia. Am J Trop Med Hyg. 1979 May;28(3):559-63. PubMed PMID: 453449.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481