

Boc-Val-Cit-PAB

Chemical Properties

CAS No. : 870487-09-5


Formula: C₂₃H₃₇N₅O₆

Molecular Weight: 479.57

Keep away from direct sunlight

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Boc-Val-Cit-PAB is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs).
Targets(IC50)	ADC Linker, PROTAC Linker
In vitro	ADCs consist of an antibody attached to an ADC cytotoxin through an ADC linker.

Solubility Information

Solubility	DMSO: 10 mM, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0852 mL	10.426 mL	20.852 mL
5 mM	0.417 mL	2.0852 mL	4.1704 mL
10 mM	0.2085 mL	1.0426 mL	2.0852 mL
50 mM	0.0417 mL	0.2085 mL	0.417 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Beck A, et al. Strategies and challenges for the next generation of antibody-drug conjugates. Nat Rev Drug Discov. 2017 May;16(5):315-337.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481