


5-(Trifluoromethyl)uridine

Chemical Properties

CAS No. :	21618-67-7
Formula:	C10H11F3N2O6
Molecular Weight:	312.2
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-(Trifluoromethyl)uridine is a Nucleoside Derivative - Fluoro-modified nucleoside, 5-Modified pyrimidine nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2031 mL	16.0154 mL	32.0307 mL
5 mM	0.6406 mL	3.2031 mL	6.4061 mL
10 mM	0.3203 mL	1.6015 mL	3.2031 mL
50 mM	0.0641 mL	0.3203 mL	0.6406 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481