


Stachybotrin E

Chemical Properties

CAS No. :	1510813-72-5
Formula:	C ₂₉ H ₃₉ N ₀ O ₈
Molecular Weight:	529.63
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Stachybotrin E is a useful organic compound for research related to life sciences. The catalog number is T125775 and the CAS number is 1510813-72-5.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8881 mL	9.4406 mL	18.8811 mL
5 mM	0.3776 mL	1.8881 mL	3.7762 mL
10 mM	0.1888 mL	0.9441 mL	1.8881 mL
50 mM	0.0378 mL	0.1888 mL	0.3776 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481