


Boditretninib

Chemical Properties

CAS No. :	1940165-80-9
Formula:	C ₂₃ H ₂₄ F ₂ N ₆ O
Molecular Weight:	438.47
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Boditretninib, a potent tyrosine kinase inhibitor, functions as an antineoplastic agent and is valuable in researching cancer, inflammation, neurodegenerative diseases, and certain infectious diseases.
Targets(IC50)	Others,c-Met/HGFR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2807 mL	11.4033 mL	22.8066 mL
5 mM	0.4561 mL	2.2807 mL	4.5613 mL
10 mM	0.2281 mL	1.1403 mL	2.2807 mL
50 mM	0.0456 mL	0.2281 mL	0.4561 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481