

Anti-Spike-RBD mAb

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Spike-RBD mAb is a CHO cell-derived human monoclonal IgG1 antibody that potentially treats SARS-CoV-2 by blocking the interaction between the spike protein and ACE2[1].
Targets(IC50)	SARS-CoV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481