

## (S)-1-(4-chlorophenyl)ethylamine

## Chemical Properties

CAS No. : 4187-56-8


Formula: C<sub>8</sub>H<sub>10</sub>ClN

Molecular Weight: 155.63

Store under nitrogen

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(S)-1-(4-chlorophenyl)ethylamine ,with CAS No. 4187-56-8, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. (S)-1-(4-chlorophenyl)ethylamine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.4255 mL	32.1275 mL	64.255 mL
5 mM	1.2851 mL	6.4255 mL	12.851 mL
10 mM	0.6425 mL	3.2127 mL	6.4255 mL
50 mM	0.1285 mL	0.6425 mL	1.2851 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481