

Compound T129958

Chemical Properties


CAS No. :

Formula: C₃₀H₃₄O₁₂

Molecular Weight: 586.59

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2S,3R,4S,5S,6R)-6-(((2R,3R,4R,5S,6S)-3,5-dihydroxy-6-methyl-4-(((2E)-3-phenylprop-2-enoyl)oxy)oxan-2-yl]oxy)methyl)-3,4,5-trihydroxyoxan-2-yl (2E)-3-phenylprop-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7048 mL	8.5238 mL	17.0477 mL
5 mM	0.341 mL	1.7048 mL	3.4095 mL
10 mM	0.1705 mL	0.8524 mL	1.7048 mL
50 mM	0.0341 mL	0.1705 mL	0.341 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481