


Disperse Fast Yellow 4K

Chemical Properties

CAS No. :	6300-37-4
Formula:	C ₁₉ H ₁₆ N ₄ O
Molecular Weight:	316.36
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Disperse Fast Yellow 4K (Disperse Yellow) is a color additive used in dyes.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 3.17 mg/mL (10.02 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.161 mL	15.8048 mL	31.6096 mL
5 mM	0.6322 mL	3.161 mL	6.3219 mL
10 mM	0.3161 mL	1.5805 mL	3.161 mL
50 mM	0.0632 mL	0.3161 mL	0.6322 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481