


Epiyangambin

Chemical Properties

CAS No. :	24192-64-1
Formula:	C ₂₄ H ₃₀ O ₈
Molecular Weight:	446.49
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Epiyangambin is a phytotoxic lignan, it shows inhibition of seed germination of <i>Agrostis stolonifera</i> cv. penncross (Poaceae) and inhibition of development of <i>Lactuca sativa</i> L. (Asteraceae) seedlings in a microassay using 24-well plates. Epiyangambin exhibits significant inhibition of the transformation of murine epidermal JB6 cells, with the IC ₅₀ value of 0.15 microg/mL. (+)-Epiyangambin also has anti-platelet aggregation activity.
Targets(IC ₅₀)	PAFR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2397 mL	11.1985 mL	22.3969 mL
5 mM	0.4479 mL	2.2397 mL	4.4794 mL
10 mM	0.224 mL	1.1198 mL	2.2397 mL
50 mM	0.0448 mL	0.224 mL	0.4479 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481