

24-Hydroxymomordicine III

Chemical Properties

CAS No. :

Formula: C₃₆H₅₈O₁₀

Molecular Weight: 650.9

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

24-Hydroxymomordicine III can be used in related research in the field of life sciences and its product number is TN6544.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5363 mL	7.6817 mL	15.3633 mL
5 mM	0.3073 mL	1.5363 mL	3.0727 mL
10 mM	0.1536 mL	0.7682 mL	1.5363 mL
50 mM	0.0307 mL	0.1536 mL	0.3073 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481