


## 4,18-Dihydro-4-hydroxysapriearine

## Chemical Properties

CAS No. :	2202760-88-9
Formula:	C <sub>20</sub> H <sub>26</sub> O <sub>3</sub>
Molecular Weight:	314.425
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4,18-Dihydro-4-hydroxysapriearine is a natural product for research related to life sciences. The catalog number is TN5591 and the CAS number is 2202760-88-9.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1804 mL	15.9018 mL	31.8036 mL
5 mM	0.6361 mL	3.1804 mL	6.3607 mL
10 mM	0.318 mL	1.5902 mL	3.1804 mL
50 mM	0.0636 mL	0.318 mL	0.6361 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481