


RWJ-60475

Chemical Properties

CAS No. : 204130-08-5
 Formula: C₁₃H₁₁BrNO₇P
 Molecular Weight: 404.11
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	RWJ-60475 is a novel CD45 phosphatase inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4746 mL	12.3729 mL	24.7457 mL
5 mM	0.4949 mL	2.4746 mL	4.9491 mL
10 mM	0.2475 mL	1.2373 mL	2.4746 mL
50 mM	0.0495 mL	0.2475 mL	0.4949 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481