

FD-838

Chemical Properties


CAS No. : 110341-78-1

Formula: C₂₂H₂₁N₀O₇

Molecular Weight: 411.4

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	FD-838 shows good antileishmanial and moderate anticancer activities, it can moderately inhibit the growth of cultured P388 and HL-60 cell lines. FD-838 has anti-fungal activity, it can significantly inhibit the growth of two plant fungal pathogens <i>Botrytis cinerea</i> and <i>Glomerella cingulata</i> with a minimum inhibitory concentration of 6.25 μ M for each, similar to that of the positive fungicide, carbendazim.
Targets(IC50)	Antifection

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4307 mL	12.1536 mL	24.3072 mL
5 mM	0.4861 mL	2.4307 mL	4.8614 mL
10 mM	0.2431 mL	1.2154 mL	2.4307 mL
50 mM	0.0486 mL	0.2431 mL	0.4861 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481