

Ald-Ph-amido-PEG3-C2-Pfp ester

Chemical Properties

CAS No. : 2101206-21-5


Formula: C₂₃H₂₂F₅N₇O₇

Molecular Weight: 519.42

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ald-Ph-amido-PEG3-C2-Pfp ester is a noncleavable antibody-drug conjugate (ADC) linker, classified as a polyethylene glycol (PEG) linker.
Targets(IC50)	ADC Linker,PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9252 mL	9.6261 mL	19.2522 mL
5 mM	0.385 mL	1.9252 mL	3.8504 mL
10 mM	0.1925 mL	0.9626 mL	1.9252 mL
50 mM	0.0385 mL	0.1925 mL	0.385 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481