


Z 6031

## Chemical Properties

CAS No. : 31681-13-7  
 Formula: C<sub>14</sub>H<sub>30</sub>ClNO<sub>5</sub>Si  
 Molecular Weight: 355.93  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Z 6031 is a biochemical.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8095 mL	14.0477 mL	28.0954 mL
5 mM	0.5619 mL	2.8095 mL	5.6191 mL
10 mM	0.281 mL	1.4048 mL	2.8095 mL
50 mM	0.0562 mL	0.281 mL	0.5619 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481