


U 34599

Chemical Properties

CAS No. : 38150-37-7
 Formula: C₁₈H₁₅ClN₄
 Molecular Weight: 322.79
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	U 34599 stimulates the production of tissue-type plasminogen activator.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.098 mL	15.4899 mL	30.9799 mL
5 mM	0.6196 mL	3.098 mL	6.196 mL
10 mM	0.3098 mL	1.549 mL	3.098 mL
50 mM	0.062 mL	0.3098 mL	0.6196 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481