

Compound T129311

Chemical Properties


CAS No. :

Formula: C₁₉H₂₈O₈

Molecular Weight: 384.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(4E,6E)-7-[3a,4,5-trihydroxy-7-(hydroxymethyl)-2-(2-hydroxypropan-2-yl)-2,3,3a,4,5,7a-hexahydro-1-benzofuran-6-yl]hepta-4,6-dienoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6013 mL	13.0063 mL	26.0125 mL
5 mM	0.5203 mL	2.6013 mL	5.2025 mL
10 mM	0.2601 mL	1.3006 mL	2.6013 mL
50 mM	0.052 mL	0.2601 mL	0.5203 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481