


Talsupram hydrochloride

Chemical Properties

CAS No. :	25487-28-9
Formula:	C ₂₀ H ₂₆ ClNS
Molecular Weight:	347.95
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	noradrenalin re-uptake inhibitor
Targets(IC50)	Others,Adrenergic Receptor

Solubility Information

Solubility	H ₂ O: <6.96 mg/mL,Sonication is recommended. DMSO: <34.8 mg/mL,Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.874 mL	14.3699 mL	28.7398 mL
5 mM	0.5748 mL	2.874 mL	5.748 mL
10 mM	0.2874 mL	1.437 mL	2.874 mL
50 mM	0.0575 mL	0.2874 mL	0.5748 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481