

4-Pyridylcarbinol

Chemical Properties


CAS No. : 586-95-8

Formula: C₆H₇NO

Molecular Weight: 109.12

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-Pyridylcarbinol, with CAS No. 586-95-8, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 4-Pyridylcarbinol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	9.1642 mL	45.8211 mL	91.6422 mL
5 mM	1.8328 mL	9.1642 mL	18.3284 mL
10 mM	0.9164 mL	4.5821 mL	9.1642 mL
50 mM	0.1833 mL	0.9164 mL	1.8328 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481