

12,13-dihydroxy-12-methyltetradecanoic acid

Chemical Properties

CAS No. :

Formula: C₁₅H₃₀O₄

Molecular Weight: 274.4

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	12,13-dihydroxy-12-methyltetradecanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6443 mL	18.2216 mL	36.4431 mL
5 mM	0.7289 mL	3.6443 mL	7.2886 mL
10 mM	0.3644 mL	1.8222 mL	3.6443 mL
50 mM	0.0729 mL	0.3644 mL	0.7289 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481