


L-Lysine, N6,N6'-(1,8-dioxo-1,8-octanediyl)bis-

Chemical Properties

CAS No. : 675587-75-4
 Formula: C₂₀H₃₈N₄O₆
 Molecular Weight: 430.54
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	L-Lysine, N6,N6'-(1,8-dioxo-1,8-octanediyl)bis- is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3227 mL	11.6133 mL	23.2266 mL
5 mM	0.4645 mL	2.3227 mL	4.6453 mL
10 mM	0.2323 mL	1.1613 mL	2.3227 mL
50 mM	0.0465 mL	0.2323 mL	0.4645 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481