

Compound N003-0005

Chemical Properties


CAS No. :

Formula: C₂₀H₃₆O₂

Molecular Weight: 308.506

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N003-0005 is a useful organic compound for research related to life sciences and the catalog number is T131383.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2414 mL	16.2069 mL	32.4139 mL
5 mM	0.6483 mL	3.2414 mL	6.4828 mL
10 mM	0.3241 mL	1.6207 mL	3.2414 mL
50 mM	0.0648 mL	0.3241 mL	0.6483 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481