

Compound T130761

Chemical Properties


CAS No. :

Formula: C₄₅H₃₆O₁₈

Molecular Weight: 864.77

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	5,13-bis(3,4-dihydroxyphenyl)-7-[2-(3,4-dihydroxyphenyl)-3,5,7-trihydroxy-3,4-dihydro-2H-1-benzopyran-8-yl]-4,12,14-trioxapentacyclo[11.7.1.0 ^{2,11} .0 ^{3,8} .0 ^{15,20}]henicosa-2,8,10,15,17,19-hexaene-6,9,17,19,21-pentol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1564 mL	5.7819 mL	11.5638 mL
5 mM	0.2313 mL	1.1564 mL	2.3128 mL
10 mM	0.1156 mL	0.5782 mL	1.1564 mL
50 mM	0.0231 mL	0.1156 mL	0.2313 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481