


Melagatran monohydrate

Chemical Properties

CAS No. : 318245-80-6
 Formula: C₂₂H₃₃N₅O₅
 Molecular Weight: 447.53
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Melagatran monohydrate is a thrombin inhibitor. Melagatran prevents tissue factor expression in human platelet-monocyte heterotypic complexes. Melagatran reduces advanced atherosclerotic lesion size and may promote plaque stability in apolipoprotein E-deficient mice.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2345 mL	11.1724 mL	22.3449 mL
5 mM	0.4469 mL	2.2345 mL	4.469 mL
10 mM	0.2234 mL	1.1172 mL	2.2345 mL
50 mM	0.0447 mL	0.2234 mL	0.4469 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481