


Cabozantinib S-malate

Chemical Properties

CAS No. :	1140909-48-3
Formula:	C32H30FN3O10
Molecular Weight:	635.59
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Cabozantinib S-malate (XL184) is the s-malate salt form of cabozantinib, an orally bioavailable, small molecule receptor tyrosine kinase (RTK) inhibitor with potential antineoplastic activity.
Targets(IC50)	Apoptosis,c-Kit,c-Met/HGFR,TAM Receptor,VEGFR
In vitro	Cabozantinib has weak inhibitory activity against RON and PDGFR β with IC50 of 124 nM and 234 nM, respectively, and has low activity against FGFR1 with IC50 of 5.294 μ M. [1] Cabozantinib at low concentration (0.1-0.5 μ M) is sufficient to induce marked inhibition of constitutive and inducible Met phosphorylation and its resultant downstream signaling in MPNST cells, and inhibit HGF-induced MPNST cell migration and invasion. Cabozantinib also induces marked inhibition of Met and VEGFR2 phosphorylation in cytokine-stimulated human umbilical vein endothelial cells (HUVECs). Although Cabozantinib has no significant effect on MPNST cell growth at 0.1 μ M, Cabozantinib at 5-10 μ M significantly inhibits the MPNST cell growth. [2]
In vivo	Cabozantinib treatment at 30 mg/kg in RIP-Tag2 mice with spontaneous pancreatic islet tumors disrupts 83% of the tumor vasculature, reduces pericytes and empty basement membrane sleeves, causes widespread intratumoral hypoxia and extensive tumor cell apoptosis, and slows regrowth of the tumor vasculature after drug withdrawal, more significantly compared with XL999 that blocks VEGFR but not c-Met, leading to only 43% reduction in vascularity, suggesting that concurrent inhibition of VEGFR and other functionally relevant receptor tyrosine kinases (RTK) amplifies angiogenesis inhibition. Cabozantinib also decreases invasiveness of primary tumors and reduces metastasis. [1] Cabozantinib at 30 mg/kg/day significantly abrogates human MPNST xenografts growth and metastasis in SCID mice. [2] Administration of Cabozantinib induces dose-dependent inhibition of tumor growth in breast, lung, and glioma tumor models, in association with decreased tumor and endothelial cell proliferation as well as increased apoptosis. A single oral dose of Cabozantinib is sufficient to induce sustained tumor growth inhibition in MDA-MB-231 tumor-bearing mice and C6 tumor-bearing rats at 100 mg/kg and 10 mg/kg, respectively. [3]
Kinase Assay	The inhibition profile of cabozantinib against a broad panel of 270 human kinases is determined using luciferase-coupled chemiluminescence, 33 P-phosphoryl transfer, or AlphaScreen technology. Recombinant human full-length, glutathione S-transferase tag, or histidine tag fusion proteins are used, and half maximal inhibitory concentration (IC50) values are determined by measuring phosphorylation of peptide substrate poly

A DRUG SCREENING EXPERT

Kinase Assay	(Glu, Tyr) at ATP concentrations at or below the K_m for each respective kinase. The mechanism of kinase inhibition is evaluated using the AlphaScreen Assay by determining the IC_{50} values over a range of ATP concentrations.
Cell Research	Cells are exposed to various concentrations of Cabozantinib for 48 hours. Cell growth is determined by MTS assays using CellTiter96 Aqueous Non-Radioactive Cell Proliferation Assay kit. Absorbance is measured at a wavelength of 490 nm, and the absorbance values of treated cells are presented as a percentage of the absorbance of untreated cells. (Only for Reference)

Solubility Information

Solubility	H2O: < 1 mg/mL (insoluble or slightly soluble), DMSO: 16.67 mg/mL (26.23 mM), Sonication is recommended. Ethanol: < 1 mg/mL (insoluble or slightly soluble), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5733 mL	7.8667 mL	15.7334 mL
5 mM	0.3147 mL	1.5733 mL	3.1467 mL
10 mM	0.1573 mL	0.7867 mL	1.5733 mL
50 mM	0.0315 mL	0.1573 mL	0.3147 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- You WK, et al. Cancer Res, 2011, 71(14), 4758-4768.
- Torres KE, et al. Clin Cancer Res, 2011, 17(12), 3943-3955.
- Yakes FM, et al. Mol Cancer Ther, 2011, 10(12), 2298-2308.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481