


## 3-Amino-5-methylisoxazole

## Chemical Properties

CAS No. :	1072-67-9
Formula:	C <sub>4</sub> H <sub>6</sub> N <sub>2</sub> O
Molecular Weight:	98.1
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Amino-5-methylisoxazole is pharmaceutical intermediates, used in the production of sulfonamide drugs.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: 55 mg/mL (560.65 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween-80+45% Saline: 5 mg/mL (50.97 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	10.1937 mL	50.9684 mL	101.9368 mL
5 mM	2.0387 mL	10.1937 mL	20.3874 mL
10 mM	1.0194 mL	5.0968 mL	10.1937 mL
50 mM	0.2039 mL	1.0194 mL	2.0387 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

M J Rieder et al. Journal of Pharmacology and Experimental Therapeutics, 244(2), 724-728 (1988)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481