

m-PEG15-NHS ester

Chemical Properties

CAS No. :

Formula: C39H73NO20

Molecular Weight: 875.99

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	m-PEG15-NHS ester is a PEG-based linker for PROTACs that joins two essential ligands, crucial for forming PROTAC molecules, and enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands joined by a linker, one binding to an E3 ubiquitin ligase and the other to the target protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1416 mL	5.7078 mL	11.4157 mL
5 mM	0.2283 mL	1.1416 mL	2.2831 mL
10 mM	0.1142 mL	0.5708 mL	1.1416 mL
50 mM	0.0228 mL	0.1142 mL	0.2283 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481