

Azido-PEG3-NHS ester

Chemical Properties

CAS No. : 1092654-47-1


Formula: C12H18N4O7

Molecular Weight: 330.29

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Azido-PEG3-NHS ester is a PEG-based linker for PROTACs [proteolysis-targeting chimeras] that joins two essential ligands, facilitating selective protein degradation via the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands connected by a linker: one binds to an E3 ubiquitin ligase and the other to a target protein. They leverage the ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0276 mL	15.1382 mL	30.2764 mL
5 mM	0.6055 mL	3.0276 mL	6.0553 mL
10 mM	0.3028 mL	1.5138 mL	3.0276 mL
50 mM	0.0606 mL	0.3028 mL	0.6055 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481