

BDM31369

Chemical Properties


CAS No. : 1001467-91-9

Formula: C₁₇H₂₃N₃O₂S

Molecular Weight: 333.45

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	BDM31369 is an activator of ethionamide.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.999 mL	14.9948 mL	29.9895 mL
5 mM	0.5998 mL	2.999 mL	5.9979 mL
10 mM	0.2999 mL	1.4995 mL	2.999 mL
50 mM	0.060 mL	0.2999 mL	0.5998 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481