


## Sarmencymarin

## Chemical Properties

CAS No. : 98633-61-5  
 Formula: C<sub>30</sub>H<sub>46</sub>O<sub>8</sub>  
 Molecular Weight: 534.69  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sarmencymarin is a cardiac glycoside.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8702 mL	9.3512 mL	18.7024 mL
5 mM	0.374 mL	1.8702 mL	3.7405 mL
10 mM	0.187 mL	0.9351 mL	1.8702 mL
50 mM	0.0374 mL	0.187 mL	0.374 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481