


XCC

Chemical Properties

CAS No. : 96865-83-7
 Formula: C19H22N4O5
 Molecular Weight: 386.4
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|--|
| Description | XCC is an Adenosine receptor antagonist. |
| Targets(IC50) | Others, Adenosine Receptor |

Solubility Information

| | |
|------------|---|
| Solubility | DMSO: <19.32 mg/mL, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble) |
|------------|---|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|------------|
| 1 mM | 2.588 mL | 12.940 mL | 25.8799 mL |
| 5 mM | 0.5176 mL | 2.588 mL | 5.176 mL |
| 10 mM | 0.2588 mL | 1.294 mL | 2.588 mL |
| 50 mM | 0.0518 mL | 0.2588 mL | 0.5176 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481