


PI3K δ / γ -IN-3

Chemical Properties

CAS No. :	2730151-31-0
Formula:	C ₂₃ H ₂₀ ClN ₉ O
Molecular Weight:	473.92
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	PI3K δ / γ -IN-3 (Compound 58) is an orally active dual inhibitor of PI3K δ (IC ₅₀ : 1 nM) and PI3K γ (IC ₅₀ : 16 nM), effective in inducing apoptosis of tumor cells and applicable to B-cell malignant tumor studies.
Targets(IC ₅₀)	Apoptosis,Others,PI3K

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1101 mL	10.5503 mL	21.1006 mL
5 mM	0.422 mL	2.1101 mL	4.2201 mL
10 mM	0.211 mL	1.055 mL	2.1101 mL
50 mM	0.0422 mL	0.211 mL	0.422 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481