

A 438079

## Chemical Properties

CAS No. : 899507-36-9  
 Formula: C<sub>13</sub>H<sub>9</sub>Cl<sub>2</sub>N<sub>5</sub>  
 Molecular Weight: 306.15  
 Storage: Store at low temperature  
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	A 438079 is a potent and selective antagonist of the P2X7 receptor (pIC <sub>50</sub> : 6.9).
Targets(IC <sub>50</sub> )	P2X Receptor
In vitro	A 438079 blocks BzATP-(10 μM) evoked changes in intracellular calcium concentrations with an IC <sub>50</sub> of 321 nM in 1321N1 cells stably expressing rat P2X7 receptors. A 438079 is also selective for the P2X7 receptor, at concentrations up to 100 μM[1].
In vivo	In neuropathic rats noxious and innocuous evoked activity of different classes of spinal neurons is reduced by A 438079 (80 μmol/kg, i.v.), and it significantly raises withdrawal thresh-olds in both the SNL and CCI models[1]. A 438079 has superior neuroprotective effects compared with an equal dose of phenobarbital (25 mg/kg)[2]. A 438079 partially but significantly prevents the 6-OHDA-induced depletion of striatal DA stores[3]. Pretreatment with A 438079 reduces nociceptive behavior scores in the HC model[4].
Kinase Assay	Human astrocytoma cells, 1321N1, are grown to stably express rat P2X7, human P2X4, P2X2a, P2X2/3, P2X1, P2Y1 and P2Y2 recombinant receptors. Agonist, BzATP, 2,3-O-(4-ben-zoylbenzoyl)-ATP or ATP-induced changes in intracellular Ca <sup>2+</sup> concentrations are assessed in all of the cell lines using the Ca <sup>2+</sup> chelating dye, Fluo-4, in conjunction with a Fluorometric Imaging Plate Reader. The cells are plated out the day before the experiment onto poly-D-lysine-coated black 96 well plates. After the agonist addition, changes in intracellular Ca <sup>2+</sup> concentrations are recorded, per second, for 3 min. Ligands are tested at 11 half-log concentrations from 0.1nM to 100 μM. BzATP or ATP concentrations corresponds to the EC <sub>70</sub> values for each receptor to enable comparison of antagonist potencies across the multiple P2 receptor subtypes. A 438079 is added to the cell plate and fluorescence data are collected for 3 min before the addition of agonist, subsequently, data are then collected for another 2 min. The pEC <sub>50</sub> or pIC <sub>50</sub> values are derived from a single curve fit [1].
Animal Research	To confirm A 438079 reach the brain after systemic administration, P10 rat pups are injected with 5 mg/kg A 438079 and killed either 10 min, 30 min, or 2 h later (n=4 per group). Blood samples are centrifuged at 1000×g for 10 min to isolate the plasma. Samples are analyzed using LC-MS/MS by a service provider. Briefly, protein is precipitated from 50 μL aliquots of the individual plasma or brain tissue homogenate, and A 438079 is quantified by LC-MS/MS from a five-point standard curve [2].

## Solubility Information

Solubility	DMSO: 145 mg/mL (473.62 mM),Sonication is recommended. H2O: 0.2 mg/mL (0.65 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2664 mL	16.3319 mL	32.6637 mL
5 mM	0.6533 mL	3.2664 mL	6.5327 mL
10 mM	0.3266 mL	1.6332 mL	3.2664 mL
50 mM	0.0653 mL	0.3266 mL	0.6533 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

- McGaraughty S, et al. P2X7-related modulation of pathological nociception in rats. *Neuroscience*. 2007 Jun 8;146(4):1817-28.
- Mesuret G, et al. *CNS Neurosci Ther*. 2014 Jun;20(6):556-64.
- Marcellino D, et al. On the role of P2X(7) receptors in dopamine nerve cell degeneration in a rat model of Parkinson's disease: studies with the P2X(7) receptor antagonist A-438079. *J Neural Transm (Vienna)*. 2010 Jun;117(6):681-7.
- Martins JP, et al. The role of P2X7 purinergic receptors in inflammatory and nociceptive changes accompanying cyclophosphamide-induced haemorrhagic cystitis in mice. *Br J Pharmacol*. 2012 Jan;165(1):183-96.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481