


HJ-PI01

## Chemical Properties

CAS No. : 6192-43-4  
 Formula: C<sub>14</sub>H<sub>11</sub>NO<sub>2</sub>  
 Molecular Weight: 225.24  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	HJ-PI01 (HJ-PI01) is a Pim-2 inhibitor. HJ-PI01 (HJ-PI01) induces apoptosis and autophagic cell death in triple-negative human breast cancer.
Targets(IC50)	Pim

## Solubility Information

Solubility	DMSO: 65 mg/mL (288.58 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4397 mL	22.1985 mL	44.3971 mL
5 mM	0.8879 mL	4.4397 mL	8.8794 mL
10 mM	0.444 mL	2.2199 mL	4.4397 mL
50 mM	0.0888 mL	0.444 mL	0.8879 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Zhao YQ, Yin YQ, Liu J, Wang GH, Huang J, Zhu LJ, Wang JH. Characterization of HJ-PI01 as a novel Pim-2 inhibitor that induces apoptosis and autophagic cell death in triple-negative human breast cancer. Acta Pharmacol Sin. 2016 Sep;37(9):1237-50.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481