

2-{2-[2-(dihydroxy-dimethyltetradecanamido)-3-methylbutanamido]-3-hydroxybutanamido} prop-enoic acid

Chemical Properties

CAS No. :

Formula: C₂₈H₅₁N₃O₈

Molecular Weight: 557.73

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	2-{2-[2-(3,5-dihydroxy-2,4-dimethyltetradecanamido)-3-methylbutanamido]-3-hydroxybutanamido}prop-2-enoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.793 mL	8.9649 mL	17.9298 mL
5 mM	0.3586 mL	1.793 mL	3.586 mL
10 mM	0.1793 mL	0.8965 mL	1.793 mL
50 mM	0.0359 mL	0.1793 mL	0.3586 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481