

## Compound T128083

## Chemical Properties

CAS No. :

Formula: C58H104O28

Molecular Weight: 1249.44

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	5-[[5-({3-hydroxy-5-[(3-hydroxy-5-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]decanoyl)oxy]decanoyl)oxy)-3-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]decanoyl]oxy)-3-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]decanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8004 mL	4.0018 mL	8.0036 mL
5 mM	0.1601 mL	0.8004 mL	1.6007 mL
10 mM	0.080 mL	0.4002 mL	0.8004 mL
50 mM	0.016 mL	0.080 mL	0.1601 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481