

WZ4141R

Chemical Properties


CAS No. :

Formula: C₁₆H₁₅N₃O₂

Molecular Weight: 281.31

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	WZ4141R is an intermediate in compound synthesis.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5548 mL	17.774 mL	35.548 mL
5 mM	0.711 mL	3.5548 mL	7.1096 mL
10 mM	0.3555 mL	1.7774 mL	3.5548 mL
50 mM	0.0711 mL	0.3555 mL	0.711 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481