


Soterenol

Chemical Properties

CAS No. :	13642-52-9
Formula:	C ₁₂ H ₂₀ N ₂ O ₄ S
Molecular Weight:	288.36
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Soterenol is an adrenergic beta-2 agonist used as bronchodilator and has a number of gastrointestinal side effects.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4679 mL	17.3394 mL	34.6789 mL
5 mM	0.6936 mL	3.4679 mL	6.9358 mL
10 mM	0.3468 mL	1.7339 mL	3.4679 mL
50 mM	0.0694 mL	0.3468 mL	0.6936 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481