


## Locicortolone dicibate

## Chemical Properties

CAS No. :	78467-68-2
Formula:	C <sub>36</sub> H <sub>50</sub> Cl <sub>2</sub> O <sub>5</sub>
Molecular Weight:	633.68
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Locicortolone dicibate is a synthetic glucocorticoid corticosteroid which was never marketed.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5781 mL	7.8904 mL	15.7808 mL
5 mM	0.3156 mL	1.5781 mL	3.1562 mL
10 mM	0.1578 mL	0.789 mL	1.5781 mL
50 mM	0.0316 mL	0.1578 mL	0.3156 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481