


SKF-83822

Chemical Properties

CAS No. : 74115-08-5
 Formula: C₂₀H₂₂ClNO₂
 Molecular Weight: 343.85
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	SKF-83822 is an agonist of dopamine D1-like receptor.
Targets(IC50)	Others,Dopamine Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9082 mL	14.5412 mL	29.0824 mL
5 mM	0.5816 mL	2.9082 mL	5.8165 mL
10 mM	0.2908 mL	1.4541 mL	2.9082 mL
50 mM	0.0582 mL	0.2908 mL	0.5816 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Paci P, Gabriele S, Ris L. A new method allowing long-term potentiation recordings in hippocampal organotypic slices. *Brain Behav.* 2017 Apr 12;7(5):e00692. doi: 10.1002/brb3.692. eCollection 2017 May. PubMed PMID: 28523233; PubMed Central PMCID: PMC5434196.

Abraham AD, Neve KA, Lattal KM. Activation of D1/5 Dopamine Receptors: A Common Mechanism for Enhancing Extinction of Fear and Reward-Seeking Behaviors. *Neuropsychopharmacology.* 2016 Jul;41(8):2072-81. doi: 10.1038/npp.2016.5. Epub 2016 Jan 14. PubMed PMID: 26763483; PubMed Central PMCID: PMC4908654.

Komal P, Estakhr J, Kamran M, Renda A, Nashmi R. cAMP-dependent protein kinase inhibits $\alpha 7$ nicotinic receptor activity in layer 1 cortical interneurons through activation of D1/D5 dopamine receptors. *J Physiol.* 2015 Aug 15; 593(16):3513-32. doi: 10.1113/JP270469. Epub 2015 Jun 25. PubMed PMID: 25990637; PubMed Central PMCID: PMC4560582.

Gildea JJ, Shah IT, Van Sciver RE, Israel JA, Enzensperger C, McGrath HE, Jose PA, Felder RA. The cooperative roles of the dopamine receptors, D1R and D5R, on the regulation of renal sodium transport. *Kidney Int.* 2014 Jul;86(1):118-26. doi: 10.1038/ki.2014.5. Epub 2014 Feb 19. PubMed PMID: 24552847; PubMed Central PMCID: PMC4077925.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481