

Boc-PEG4-C2-NHS ester

Chemical Properties

CAS No. : 1835759-75-5


Formula: C₂₀H₃₃N₁O₁₀

Molecular Weight: 447.48

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Boc-PEG4-C2-NHS ester is a PEG-based linker for PROTACs that joins two essential ligands to form PROTAC molecules, enabling selective protein degradation via the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs, comprising two distinct ligands connected by a linker, utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins; one ligand binds to an E3 ubiquitin ligase while the other binds to the target protein [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2347 mL	11.1737 mL	22.3474 mL
5 mM	0.4469 mL	2.2347 mL	4.4695 mL
10 mM	0.2235 mL	1.1174 mL	2.2347 mL
50 mM	0.0447 mL	0.2235 mL	0.4469 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481