


## Epostatin

## Chemical Properties

CAS No. :	181372-99-6
Formula:	C <sub>23</sub> H <sub>33</sub> N <sub>3</sub> O <sub>5</sub>
Molecular Weight:	431.53
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Epostatin is a new dipeptidyl peptidase II (DPP-II, EC 3.4.14.2) inhibitor.
Targets(IC50)	Others,Proteasome

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3173 mL	11.5867 mL	23.1734 mL
5 mM	0.4635 mL	2.3173 mL	4.6347 mL
10 mM	0.2317 mL	1.1587 mL	2.3173 mL
50 mM	0.0463 mL	0.2317 mL	0.4635 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

### Reference

- Grzegorzewska AE, Winnicka H, Warchoń W, Mostowska A, Jagodziński PP. Correlations of indoleamine 2,3-dioxygenase, interferon- $\lambda$ 3, and anti-HBs antibodies in hemodialysis patients. *Vaccine*. 2018 Jun 20. pii: S0264-410X(18)30862-4. doi: 10.1016/j.vaccine.2018.06.034. [Epub ahead of print] PubMed PMID: 29935858.
- Meijssen JJ, Rammos A, Campbell A, Hayward C, Porteous DJ, Deary IJ, Marioni RE, Nicodemus KK. Using tree-based methods for detection of gene-gene interactions in the presence of a polygenic signal: simulation study with application to educational attainment in the Generation Scotland Cohort Study. *Bioinformatics*. 2018 Jun 19. doi: 10.1093/bioinformatics/bty462. [Epub ahead of print] PubMed PMID: 29931044.
- Chen B, Chen J, Du Q, Zhou D, Wang L, Xie J, Li Y, Zhang D. Genetic variants in microRNA biogenesis genes as novel indicators for secondary growth in *Populus*. *New Phytol*. 2018 Jun 19. doi: 10.1111/nph.15262. [Epub ahead of print] PubMed PMID: 29916214.
- Yu Y, Assmann SM. Inter-relationships between the heterotrimeric G $\beta$  subunit AGB1, the RLK FERONIA and RALF1 in salinity response. *Plant Cell Environ*. 2018 Jun 15. doi: 10.1111/pce.13370. [Epub ahead of print] PubMed PMID: 29907954.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481