

Compound N032-0003

Chemical Properties


CAS No. :

Formula: C₃₂H₅₂O₃

Molecular Weight: 484.765

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N032-0003 is a useful organic compound for research related to life sciences and the catalog number is T131300.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0628 mL	10.3142 mL	20.6283 mL
5 mM	0.4126 mL	2.0628 mL	4.1257 mL
10 mM	0.2063 mL	1.0314 mL	2.0628 mL
50 mM	0.0413 mL	0.2063 mL	0.4126 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481