

Compound T130310

Chemical Properties


CAS No. :

Formula: C27H30O13

Molecular Weight: 562.52

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

10-(4-hydroxy-3,5-dimethoxyphenyl)-16-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]-4,6,13-trioxatetracyclo[7.7.0.0^{3,7}.0^{11,15}]hexadeca-1,3(7),8-trien-12-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7777 mL	8.8886 mL	17.7771 mL
5 mM	0.3555 mL	1.7777 mL	3.5554 mL
10 mM	0.1778 mL	0.8889 mL	1.7777 mL
50 mM	0.0356 mL	0.1778 mL	0.3555 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481