


cIAP1 Ligand-Linker Conjugates 15 hydrochloride

Chemical Properties

CAS No. :	1225383-36-7
Formula:	C37H48ClN3O8
Molecular Weight:	698.25
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	cIAP1 Ligand-Linker Conjugates 15 hydrochloride is a chemical compound that combines an IAP ligand for the E3 ubiquitin ligase and a PROTAC linker, useful in the development of SNIPERs [1].
Targets(IC50)	Others,E3 Ligase Ligand-Linker Conjugates,PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4322 mL	7.1608 mL	14.3215 mL
5 mM	0.2864 mL	1.4322 mL	2.8643 mL
10 mM	0.1432 mL	0.7161 mL	1.4322 mL
50 mM	0.0286 mL	0.1432 mL	0.2864 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Itoh Y, et al. Protein knockdown using methyl bestatin-ligand hybrid molecules: design and synthesis of inducers of ubiquitination-mediated degradation of cellular retinoic acid-binding proteins. J Am Chem Soc. 2010 Apr 28;132(16):5820-6.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481